
The CrossFix® II System


All-Suture, All-Inside, Meniscal Repair



Fast, Easy Reproducible Technique Stronger, Stiffer Needles** Improved Suture Strength** Instantaneous Mattress Stitch



All inside Technique. Inside out Results.



Suture only device eliminates concerns associated with hard implants

Curved and straight delivery needles allow for access to multiple tear sites

High strength suture for easy knot sliding and high knot breakage strength

Sharper, stiffer needles** improves delivery and consistency

Incorporated depth limiter provides for needle penetration control

CrossFix II Meniscal Repair System

Description	Part Number
CrossFix II Disposable Procedure Kit (Probe, Introducer and Knot Pusher/Suture Cutter)	CM-8000
CrossFix II Meniscal Repair Device, Straight	CM-8001
CrossFix II Meniscal Repair Device, Up Curve	CM-8002



CrossFix II Disposable Procedure Kit

- *Barber, FA, et al. "Biomechanical Testing of New Meniscal Repair Techniques Containing Ultra High Molecular Weight Polyethylene Suture." Arthroscopy. 2009; Vol. 25, No 9: pp. 959-967.
- **Compared to the original CrossFix iteration; bench testing is not necessarily indicative of clinical performance

This material is intended for health care professionals. Distribution to any other recipient is prohibited. All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated. This material must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet. Check for country product clearances and reference product specific instructions for use. For complete product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see the package insert and Zimmer Biomet's website. Not for distribution in France.

©2018 Zimmer Biomet





Legal Manufacturer Cayenne Medical Inc. 16597 N. 92nd St Suite 101 Scotsdale, AZ 85260 USA

www.zimmerbiomet.com